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Tel.: +49(0)7151 975099-0
Fax: +49(0)7151 975099-30

Internet: www.ekastu.de
eMail: info@ekastu.de

Information supplied by the Manufacturer for:

EKASTU-Eye Wash Bottle with funnel, 600 ml FD *	(Order no.: 177 959)
EKASTU-Eye Wash Bottle ECO with funnel, 600 ml FD **	(Order no.: 177 111)
EKASTU-Eye Wash Bottle MINI with funnel, 200 ml FD *	(Order no.: 177 259)
EKASTU-Eye Wash Bottle MINI-ECO with funnel, 200 ml FD **	(Order no.: 177 112)
EKASTU-Eye Wash Bottle ADR200, FD	(Order no.: 177 312)

DIN EN 15154-4

These Medical Devices of class I are in compliance with Regulation (EU) 2017/745, especially all essential requirements of Annex I and standard DIN EN 15154-4:2009-07. The EU Declaration of Conformity is available for download from www.ekastu.de.

Purpose

Eye wash bottles are used to rinse foreign matters out of the eye, as part of first aid. Eye wash bottles allow a continued flushing on the way to medical care also. Contraindications and side effects are not known.

General Notes

Use of an eye wash bottle presupposes knowledge of and adherence to this information brochure. Additionally, the information noted directly on the eye wash bottles must be considered. The eye wash bottles are designed exclusively for the use named in this information brochure and on the bottles. EKASTU Safety GmbH assumes liability within of the General Terms and Conditions. Insofar you do not have them, they will be send to you on request. The manufacturer is not liable for damage caused by ignoring the information. Serious incidents related to the product must be reported to the manufacturer.

Use

How to use these eye wash bottles is described on their labels step to step.

All eye wash bottles are in accordance to DIN EN 15154-4.

The first seconds after an accident, where the eyes are exposed to dangerous substances, are decisive in keeping the injury as low as possible.

Non plumbed-in eye wash units (eye wash bottles) are designed and intended to be kept in direct proximity to those people working in potentially dangerous areas. The principal duty of these items of equipment is to provide first aid.

Eye wash bottles also enable continual rinsing whilst on the way to get medical treatment.

The bottles are made of transparent polyethylene and are leak proof.

The screwed-on eye funnel (*) (**), which fits around the eye socket, guarantees safe protective contact for the eyeball.

Contaminated water drains through a separate discharge point (*) (**).

Storage and Shelf Life

Please consider the pictograms regarding storage, noted on the bottles.

EKASTU-Eye Wash Bottles shall be protected from dust and irradiation and oxidant agents.

The expiry date for correct stored and originally closed eye wash bottles is indicated on the EKASTU-Eye Wash Bottles. Do not use open or damaged eye wash bottles. Failure to do so may result in contamination and resulting infection.

Recommendation for Storage

For optimal storage of EKASTU-Eye Wash Bottles we recommend the Eye Wash Station (EPS) empty. Order no.: 277 955

Pictogram-Definition

For single use only

See information supplied by the manufacturer

Replaces the word „order number“


Medical device

Batch number

Unique Device Identification

Swiss Representative

End of shelf life

Do not use the content if packaging (bottle) is damaged

Storage temperature range

Protect against solar irradiation

Manufacturer

CE-Marking

Pharmaceutical registration number

Pictogram for eye wash device

Use the eye wash bottle in lying position

Use the eye wash bottle in upright position

Serious incidents related to the product must be reported to the manufacturer

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